

# Q-MACS Trace compact

## technical specification

The Q-MACS Trace compact is a newly developed compact, easy portable system for high sensitive trace gas measurements using a 36 m optical long path cell. It is based on the Q-MACS Basic and uses infrared absorption spectroscopy to measure absolute molecular concentrations.



### general

description	single path infrared spectrometer with IR-light source and 36 m long path cell
sensitivity	down to ppb range [1]
time resolution	down to milliseconds
size	520 mm x 600 mm x 370 mm
weight	40 kg
system version	Q-MACS Trace 3.6 - 36 m Long Path Cell

### components

parts	<ul style="list-style-type: none"> <li>▪ optical board</li> <li>▪ laptop with PCMCIA-card</li> <li>▪ 36 m Herriot cell</li> <li>▪ vacuum pump system</li> <li>▪ water cooling system</li> </ul>
-------	---

### parameter

power	<ul style="list-style-type: none"> <li>▪ 230 V, max. 2 A (switch-on current 6 A)</li> <li>▪ 115 V, max. 4 A (switch-on current 12 A)</li> </ul>
working range	+5 °C to +40 °C

### QCL

tuning method	<ul style="list-style-type: none"> <li>▪ inter pulse mode (laser sweep mode)</li> <li>▪ intra pulse mode (single pulse mode)</li> </ul>
pulse width	8 ns* ... 256 ns** * depends on the QCL and QCL-voltage used ** longer pulses on request
pulse frequency	up to 1 MHz
QCL temperature range	-35 °C to +40 °C
QCL	tested and installed

[1] depends on species, temperature and pressure

